


GAMRY REFERENCE 600: POTENTIOSTAT/GALVANOSTAT/ZRA

Overview

The Reference 600 is a high-performance, research-grade potentiostat/galvanostat/ZRA designed for fast, low-current measurements. It does well for a variety of applications such as physical electrochemistry (especially at microelectrodes), fast cyclic voltammetry, electrochemical corrosion, electrochemical noise measurements, paints and coatings, and sensors. It has a number of auxiliary input and outputs designed to help you interface or control ancillary equipment such as a rotating disc electrode. It also has a thermocouple for temperature measurements.

Electrochemical Impedance Spectroscopy

The Reference 600 comes fully equipped to perform electrochemical impedance spectroscopy. The Accuracy Contour Plot shown below provides a detailed look at the performance you can expect from your instrument in real-world situations. The results include the cell cable.



Below are additional details regarding the capabilities of the Reference 600 potentiostat. Each bullet point contains a list of the type of techniques available for the instrument to run.

- **Physical Electrochemistry** - Techniques such as cyclic voltammetry, chronoamperometry, and chronopotentiometry and derivatives of these techniques.
- **Pulse Voltammetry** - Techniques such as pulse voltammetry, square wave voltammetry, and associated stripping techniques such as anodic stripping voltammetry.
- **DC Corrosion** - Run standard DC corrosion tests such as polarization resistance, potentiodynamic, cyclic polarization, and galvanic corrosion in addition to a number of others.

- **Electrochemical Energy** - Test single-cells and stacks of various batteries, fuel cells or supercapacitors. Includes charge, discharge, cyclic charge-discharge techniques, potentiostatic, galvanostatic, self-discharge, leakage rate, and read cell voltage.
- **Electrochemical Signal Analyzer** - Designed specifically for the acquisition and analysis of time-dependent electrochemical noise signals. Cell voltage and current are continuously monitored at rates from 0.1 Hz to 1 kHz. A full featured set of analysis tools provides powerful analysis features such as statistical analysis, detrending, impedance spectra, and histogram analysis.
- **Electrochemical Frequency Modulation** - A non-destructive corrosion rate measurement technique. It allows for measurement of the corrosion rate without prior knowledge of the Tafel constants. In addition, the technique determines the Tafel constants and provides 2 internal validity checks.
- **Critical Pitting Temperature** - controls a Gamry Potentiostat, TDC4 Temperature Controller, and associated accessories to automatically measure the Critical Pitting Temperature of a material.
- **Electrochemical Noise** - A more general form of electrochemical noise testing. It is also an ECM8 Multiplexer compatible electrochemical noise software package.
- **Electrochemical Impedance Spectroscopy** - includes experimental scripts for potentiostatic, galvanostatic and hybrid impedance spectroscopy experiments in addition to single frequency techniques like Mott-Schottky. We also have our unique power-leveling multisine technique that improves signal-to-noise across the spectrum. On the analysis side, it provides tools for fitting spectra to equivalent circuit models, Kramers-Kronig transform for data validation and a graphical model editor. Our software even includes a script for EIS simulation.
- **eChemAC** - Includes full capabilities of eChemDC Toolkit plus allows electrochemical impedance spectroscopy (EIS) and EFM experiments.

Other options

- 2, 3, and 4 electrode measurements
- Electrical Isolation
 - Floating instrument: use with autoclaves, mechanical stress apparatus, or pipeline probes.
- Portable
 - Size of a chemistry textbook, weighing only 3 kgs (6.6 lbs). Easy USB 2.0 connection to a Windows computer.
- Built-In EIS
 - On-board DDS to perform EIS from 10 μ Hz to 1 MHz.
- DSP (Digital Signal Processing) Mode
 - Oversamples for improved signal-to-noise and accurate capacitance measurements.
- Current Interrupt and Positive Feedback iR Compensation
 - Gamry potentiostats and their controlling software use control loop algorithms to accurately measure and correct for uncompensated resistance.
- Auxiliary I/O
 - Control additional equipment via additional I/O interfaces: external signal input, analog voltage output, analog current output, auxiliary A/D input, and digital I/O connector.
- Warranty
 - Protected by 2-year factory service warranty.

SPECIFICATIONS:

	REFERENCE 600
Potentiostat	YES
Galvanostat	YES
Zero Resistance Ammeter	YES
Cell Connections	2, 3, or 4
Floating (Isolated from Earth)	YES
<i>WEIGHT</i>	3 kg
<i>DIMENSIONS</i>	9 (W) x 19 (H) x 27 (D) cm
<i>SYSTEM</i>	
Max Current	± 600 mA
Current Ranges	11 (60 pA - 600 mA)
Current Ranges (w/Internal Gain applied)	13 (600 fA - 600 mA)
Min Voltage Resolution	1 μV
Min Current Resolution	20 aA
Max Applied Potential	± 11 V
Rise Time	<250 ns
Noise and Ripple	<10 μV rms

Noise and Ripple (typical)	<2 μV rms
Min Time Base	3.333 μs
Max Time Base	715 s
Min Potential Step	12.5 μV
<i>EIS MEASUREMENT</i>	
Frequency Range	10 μHz - 1 MHz
EIS Accuracy	See Accuracy Contour Plot
Max AC Amplitude	3 V max 600 mA max
<i>CONTROL AMP</i>	
Compliance Voltage	± 22 V
Output Current	$> \pm 600$ mA
Speed Settings	5
Unity Gain Bandwidth	980, 260, 40, 4, 0.4 kHz
<i>ELECTROMETER</i>	
Input Impedance	$>10^{14}$ Ω
Input Current	<5 pA
Input Current (typical)	<2 pA

Bandwidth (-3dB) (typical)	>15 MHz
Common Mode Rejection Ratio	>80 dB (3 Hz) >60 dB (1 MHz)
<i>APPLIED POTENTIAL</i>	
Accuracy	$\pm 1 \text{ mV} \pm 0.2\%$ of setting
Accuracy (typical)	$\pm 375 \text{ }\mu\text{V} \pm 0.04\%$ of reading
Resolution	12.5 μV , 50 μV , 200 $\mu\text{V/bit}$
Drift	<20 $\mu\text{V}/^\circ\text{C}$
Potential Scan Range	$\pm 0.4 \text{ V}$, $\pm 1.6 \text{ V}$, $\pm 6.4 \text{ V}$
<i>MEASURED POTENTIAL</i>	
Accuracy	$\pm 1 \text{ mV} \pm 0.2\%$ of reading
Accuracy (typical)	$\pm 250 \text{ }\mu\text{V} \pm 0.05\%$ of reading
Full-Scale Ranges	12 V, 3 V, 300 mV, 30 mV
Resolution	400 μV , 100 μV , 10 μV , 1 $\mu\text{V/bit}$
Offset Range	$\pm 10 \text{ V}$
<i>CURRENT</i>	
Applied/Measured Accuracy	$\pm 10 \text{ pA} \pm 0.05\%$ of range $\pm 0.2\%$ of value (600 mA-6 nA) or 0.75% of value (600 pA) or 1.5% of value (60 pA)

Applied/Measured Resolution	0.0033% full-scale/bit
Bandwidth (-3 dB) NOTE: Bandwidth is current-range dependent.	> 10 MHz (600 mA – 600 μ A) > 1.5 MHz (60 μ A) > 0.15 MHz (6 μ A)
Stability Settings	4
Post Offset Gain	1, 10, 100
Offset Range	\pm 1X full-scale
<i>iR COMPENSATION</i>	
Mode	Current Interrupt and Positive Feedback
Min Interrupt Time	33 μ s
Max Interrupt Time	715 s
<i>AUXILIARY A/D INPUT</i>	
Range	\pm 3 V
Resolution	0.1 mV
Input Impedance	>100 k Ω or >10 G Ω
<i>AUXILIARY D/A OUTPUT</i>	
Range	0-4 V
Resolution	1 mV